



# KERN SET-06 EasyTouch Classification



1 Display of weighing class




2 Predefined/new tolerance classes

## ET Classification – Classification function


### Features

- Note: the required basis is **SET-01 ET OS** (basic package)
- The job of the classification function is to group weighing objects into specific weight classes. These weight classes can be freely determined. This function is for simple, convenient sorting of objects and avoids the user having to read out and interpret the weighing result
- **1** Display of the **classification result**: In the result, the right weight class with regard to the current weighing object will be identified in clear text. In addition to this, the current weight class will be highlighted as a coloured field and the weighing result will be shown within the weight class field, similar to a bar graph display, with a small arrow as a visual support
- **2** It is easy to enter **the tolerance classes** with individual class names (which are also automatically incremented) with lower and upper limits, possibly with a price and individual comment for each class. As an alternative, predefined class names can also be selected and given individual limits
- **Endless classes** possible: Class limits from 0 to a target value or from a target value to infinity can be selected as class limits. You can also set class limits with gaps. These gaps are treated as “class-less”


- Central **master data memory**: Classification objects can be stored in the system memory with classification profile, class name, lower and upper limits, possibly price and comment for each class, as well as ID number and ID name, image etc. By doing this, these values do not have to be constantly re-entered, but can easily be recalled from the memory. In the master data memory, you can store a possible tare value for typical packaging, box or container, which would typically be used for the item and which is then automatically subtracted from the weighing result.
- **3 ID security**: Offers the possibility of saving every weighed and stored weighing result with a unique ID number (Dynamic Object ID) and an ID name (Dynamic Object Name). The saving process can occur on a semi-automatic or fully-automatic basis and certainly every time the load is taken off the balance and then load is applied again. This means that the user does not have to press any buttons for mass storage and can work efficiently **4**
- **5 PC print function and barcode scanning function**: By operating the KERN EasyTouch App in a Windows® or Android™ environment you can use the full PC/tablet accessory infrastructure. In particular, standard Windows printers and PC label printers can print out extensive counting slips or compact adhesive labels with the count result to suit your requirements



3 Input of tolerance classes




4 ID security




5 Mass storage


# KERN SET-06 EasyTouch Classification



Print protocol - classification



Master data memory



Sound output

## Options

- 6 Save-Server** central data memory function for additional storage of all measuring data in a central, local server directory. This is where measuring data is stored from all weighing systems connected using KERN EasyTouch, as well as all installed KERN EasyTouch functions. The advantage of this, particularly for users with several weighing systems, is having all weighing data consolidated in just one database and only having to search for individual measuring data from different balances in one table. Save-Server data storage is also tamper-proof and cannot be changed, KERN SET-10
- 7 SET-261 Voice Output:** With this option, individual voice or sound files can be stored in the system for specific events. As soon as the event occurs, the system plays back the individual sound file instead of the standard sound file. In this way, for example, clear speech output can be stored for workshops for the blind, such as “Class 1”, “Type B” or “too heavy”, KERN SET-261

- Save-Cloud:** Has the same central data memory function as the Save-Server for all weighing systems connected to KERN EasyTouch. The difference is that in Save-Data Cloud, the storage location is a KERN Server, which can be accessed over the internet, instead of a server in a local network. Setting up Save-Data Cloud functionality takes place automatically and does not require a network administrator on the user side, KERN SET-10

## Technical data

- Licence model:** A license can be operated on up to four terminal devices (PCs, laptops, tablets) at the same time and independently
- User:** An unlimited number of users can be created in one license
- Balances:** You can create and operate as many balances in one licence as you want
- Communication between balance/terminal device:** Balances can communicate with the PC, laptop or tablet by serial connection, USB, Bluetooth, Ethernet or WIFI

