

Stereomicroscope Modular System - Universal stands KERN OZB-UP



Individuality, variety and flexible working through our modular construction system ► PREMIUM universal stands

Features

- With our universal stands and basic stands, as well as microscope heads and external illumination units, you can configure your microscope to your own specifications and adapt it to your application
- Thanks to the versatile, adjustable universal stands it is possible to work in the very best way in all areas with the most varied of samples
- Large universal stands are available as stand base variants as well as with the option of a clamp for the edge or the centre of a bench. Depending on the model, you have the choice of a telescopic arm stand, a jointed arm stand or a telescopic double arm universal stand with ball bearings

Technical data

• Column height: 515 mm

OZB-A5201/OZB-A5211/OZB-A5221

• Length telescopic arm: 614 mm

OZB-A5212/OZB-A5222

• Length jointed arm: 553 mm

OZB-A5203/OZB-A5213/OZB-A5223

• Length double arm: 545 mm

Model	Description
KERN	
OZB-A5201	Telescopic arm - Plate - excl. holder
OZB-A5211	Telescopic arm - Clamp Edge of bench (Range: max. 62 mm) - excl. holder
OZB-A5221	Telescopic arm – Clamp Centre of bench (hole required) – excl. holder
OZB-A5212	Jointed arm - Clamp Edge of bench (Range: max. 62 mm) - excl. holder
OZB-A5222	Jointed arm – Clamp Centre of bench (hole required) – excl. holder
OZB-A5203	Telescopic double arm with ball bearings – Plate – excl. holder
OZB-A5213	Telescopic double arm with ball bearings – Clamp Edge of bench (Range: max. 62 mm) – excl. holder
OZB-A5223	Telescopic double arm with ball bearings - Clamp Centre of bench (hole required) - excl. holder

MICROSCOPES & REFRACTOMETERS 2024

KERN Pictograms





360° rotatable microscope head



Monocular MicroscopeFor the inspection with one eye



Binocular MicroscopeFor the inspection with both eyes



Trinocular MicroscopeFor the inspection with both eyes and the additional option for the connection of a camera



Abbe Condenser

With high numerical aperture for the concentration and the focusing of light



Halogen illumination For pictures bright and rich in contrast



LED illumination

Cold, energy-saving and especially long-life illumination



Incident illuminationFor non-transparent objects



Transmitting illuminationFor transparent objects



Fluorescence illumination For stereomicroscopes



Fluorescence illumination for compound microscopes

With 100 W mercury lamp and filter



Fluorescence illumination for compound microscopes

With 3 W LED illumination and filter



Phase contrast unit

For a higher contrast



Darkfield condenser/ unit

For a higher contrast due to indirect illumination



Polarising unit
To polarise the light



Infinity system Infinity corrected optical system



Zoom magnification For stereomicroscopes



·



Auto-focus

For automatic control of the focus level



Parallel optical system For stereomicroscopes, enables fatigue-proof working



Integrated scale

In the eyepiece



SD cardFor data storage



USB 2.0 interfaceFor data transmission



USB 3.0 interface For data transmission



WIFI data interface:

For transmitting of the picture to a mobile display device



HDMI digital camera

For direct transmitting of the picture to a display device



PC software

To transfer the measurementsfrom the device to a PC.



Automatic temperature compesation

For measurements between 10 °C and 30 °C



Protection against dust and water splashes IPxx:

The type of protection is shown in the pictogram cf. DIN EN 60529:2000-09, IEC 60529:1989+A1:1999 +A2:2013



Battery operation

Ready for battery operation. The battery type is specified for each device.



Battery operation rechargeable

Prepared for a rechargeable battery operation



Plug-in power supply

230V/50Hz in standard version for EU.
On request GB, AUS or USA version.



Integrated power supply unit

Integrated in microscope. 230V/50Hz standard EU. More standards e.g. GB, AUS or USA on request.



Package shipment

The time required to manufacture the product internally is shown in days in the pictogram.



Pallet shipment

The time required to manufacture the product internally is shown in days in the pictogram.

Abbreviations

C-Mount Adapter for the connection of a

camera to a trinocular microscope

FPS Frames per second

H(S)WF High (Super) Wide Field (Eyepiece with high eye

point for wearers of glasses)

LWD Long Working Distance

N.A. Numerical Aperture

SLR camera Single-Lens Reflex camera

SWF Super Wide Field (Field number at least Ø 23 mm

for 10× eyepiece)

W.D. Working Distance

WF Wide Field (Field number up to Ø 22 mm

for 10× eyepiece)