


Easy  
Touch



# KERN EASYTOUCH


WE MAKE YOUR WEIGHING DATA DIGITAL



# WHY KERN EASYTOUCH

INDIVIDUAL. PRODUCTIVE. EFFECTIVE.

Simple, self-explanatory and intuitive operation via a clear graphical interface with user administration and individually controllable user rights. Local and centralised data storage is possible, as well as flexible connectivity to many types of balances and model series.



Do you want to store weighing data easily and securely? Do you want to support your users in their weighing processes? Do you want to reduce or prevent measuring errors and process failures? Then expand your balances with the function you need from the KERN EasyTouch software. Whether it's for

complex recipes, accurate piece counting or rapid, effective portioning – KERN EasyTouch always has the right function. Faster working, more convenient record keeping and higher levels of process security are benefits for you with KERN EasyTouch – the easy way.

# THE FEATURES AT A GLANCE

WEIGHING, COUNTING AND DATA PROCESSING.



Create master data and container data or import, save and manage from ERP.


Export master, measurement and product data, e.g. to ERP, PPS or patient file systems, or print out individually.



Recipe function with user, component and recipe management, target dosing, multiplier function and automatic correction function in the event of overdosing.



Security in sensitive and critical processes.



Data determination and checks in accordance with the Prepackaging Ordinance (Fertigpackungsverordnung, FPVO) with weighing and tolerance mode and customisable print designs for labelling.



Compliant with 21 CFR part 11 for pharmaceutical production

Tolerance weighing with graphical display, numerical or percentage tolerance and integrated master data management.



Counting objects via reference weight from master data, target counting, or piece counting/ picking individual components in parts lists.

# INDIVIDUAL

OUR DEFINITION OF EASY




For us, IoT ready means that these balances can be easily and flexibly integrated into the existing and future network of your company. Data transfer from the balance to your delivery note or as a process in your ERP system: This is the future of weighing with EasyTouch.

Error-free, efficient and with less manpower.

## A powerful weighing system is that simple, with KERN EasyTouch


WEIGHING/MEASURING DEVICE



Easy Touch by KERN



POSSIBILITIES OF PROCESSING



### YOUR ADVANTAGES

- Speed up your processes – access data in seconds, eliminating the need for manual and time-consuming data maintenance
- Avoid errors by transferring weighing data automatically
- Save yourself additional costs, as no extra software is required
- Back up your data automatically via alibi memory (factory option)\*
- Benefit from maximum flexibility and customisability at interfaces such as RS 232, USB, Bluetooth or Ethernet – easy to retrofit without recalibration\*

\* if the model series is compatible



Detailed information and more products:  
[lp.kern-sohn.com/en/easytouch](http://lp.kern-sohn.com/en/easytouch)



**Do you have questions about our KERN EasyTouch?**

I will be happy to advise you.

Wolfgang Schäfer  
Software Sales and  
Corporate Development  
Tel. +49 7433 9933-244  
[wolfgang.schaefer@kern-sohn.com](mailto:wolfgang.schaefer@kern-sohn.com)