

KERN 303-03

KERN

Weight set, E1, 1 mg - 100 g, Knob, Stainless steel polished, in wooden case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E1
Product family	303

Measuring System

Test weight - nominal value	1 mg - 100 g
Weight OIML class	E1

Construction

Shape weight	Knob
Dimension (WxDxH)	152x152x66 mm
Number of weights	21
Case of weight set	in wooden case
Material	Stainless steel polished

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (WxDxH)	152x152x66 mm
Net weight	0,9 kg
Shipping method	Parcel service
Net weight approx.	0,90 kg
Gross weight approx.	0,90 kg
Shipping weight	0,9 kg

Services (optional)

Article number for DAkkS calibration	963-202
--------------------------------------	---------