

KERN 303-044

KERN

Weight set, E1, 1 mg - 200 g, Knob, Stainless steel polished, in plastic case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E1
Product family	303

Measuring System

Test weight - nominal value	1 mg - 200 g
Weight OIML class	E1

Construction

Shape weight	Knob
Dimension (W×D×H)	250×210×75 mm
Number of weights	23
Case of weight set	in plastic case
Material	Stainless steel polished

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	250×210×75 mm
Net weight	1,27 kg
Shipping method	Parcel service
Net weight approx.	1,4 kg
Gross weight approx.	1,4 kg
Shipping weight	1,27 kg

Services (optional)

Article number for DAkkS calibration	963-203
--------------------------------------	---------