

# KERN 303-064

**KERN**

Weight set, E1, 1 mg - 1 kg, Knob, Stainless steel polished, in plastic case



## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML E1
Product family	303

## Measuring System

Test weight - nominal value	1 mg - 1 kg
Weight OIML class	E1

## Construction

Shape weight	Knob
Dimension (W×D×H)	350×320×180 mm
Number of weights	25
Case of weight set	in plastic case
Material	Stainless steel polished

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	350×320×180 mm
Net weight	3,7 kg
Shipping method	Parcel service
Net weight approx.	3,8 kg
Gross weight approx.	3,8 kg
Shipping weight	4 kg

## Services (optional)

Article number for DAkkS calibration	963-205
--------------------------------------	---------