

KERN 303-076

KERN

Weight set, E1, 1 mg - 2 kg, Knob, Stainless steel polished, in aluminium case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E1
Product family	303

Measuring System

Test weight - nominal value	1 mg - 2 kg
Weight OIML class	E1

Construction

Shape weight	Knob
Dimension (W×D×H)	278×180×138 mm
Case of weight set	in aluminium case
Material	Stainless steel polished

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	278×180×138 mm
Net weight	5,15 kg
Shipping method	Parcel service
Net weight approx.	6 kg
Gross weight approx.	6 kg
Shipping weight	5,55 kg

Services (optional)

Article number for DAkkS calibration	963-206
--------------------------------------	---------