

# KERN 303-08

**KERN**

Weight set, E1, 1 mg - 5 kg, Knob, Stainless steel polished, in wooden case



## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML E1
Product family	303

## Measuring System

Test weight - nominal value	1 mg - 5 kg
OIML-compliant	✓
Weight OIML class	E1

## Construction

Shape weight	Knob
Dimension (W×D×H)	420×200×200 mm
Number of weights	28
Case of weight set	in wooden case
Material	Stainless steel polished

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	420×200×200 mm
Net weight	17,81 kg
Shipping method	Parcel service
Net weight approx.	18 kg
Gross weight approx.	18 kg
Shipping weight	17,81 kg

## Services (optional)

Article number for DAKkS calibration	963-207
--------------------------------------	---------