

# KERN 303-09

# KERN

Weight set, E1, 1 mg - 10 kg, Knob, Stainless steel polished, in wooden case



## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML E1
Product family	303

## Measuring System

Test weight - nominal value	1 mg - 10 kg
Weight OIML class	E1

## Construction

Shape weight	Knob
Dimension (W×D×H)	540×200×240 mm
Number of weights	29
Case of weight set	in wooden case
Material	Stainless steel polished

## Packing & Shipping

Readability force [d] (N)	2 d
Dimensions packaging (W×D×H)	540×200×240 mm
Net weight	30,61 kg
Shipping method	Freight
Net weight approx.	32 kg
Gross weight approx.	34 kg
Shipping weight	32,5 kg

## Services (optional)

Article number for DAkkS calibration	963-208
--------------------------------------	---------