

KERN 304-034

KERN

Weight set, E1, 1 g - 100 g, Knob, Stainless steel polished, in plastic case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E1
Product family	304

Measuring System

Test weight - nominal value	1 g - 100 g
OIML-compliant	✓
Weight OIML class	E1

Construction

Shape weight	Knob
Dimension (W×D×H)	250×210×75 mm
Number of weights	9
Case of weight set	in plastic case
Material	Stainless steel polished

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	250×210×75 mm
Net weight	0,9 kg
Shipping method	Parcel service
Net weight approx.	0,90 kg
Gross weight approx.	0,90 kg
Shipping weight	0,9 kg

Services (optional)

Article number for DAKkS calibration	963-216
--------------------------------------	---------