

# KERN 304-034

# KERN

Weight set, E1, 1 g - 100 g, Knob, Stainless steel polished, in plastic case



## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML E1
Product family	304

## Measuring System

Test weight - nominal value	1 g - 100 g
Weight OIML class	E1

## Construction

Shape weight	Knob
Dimension (WxDxH)	250x210x75 mm
Number of weights	9
Case of weight set	in plastic case
Material	Stainless steel polished

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (WxDxH)	250x210x75 mm
Net weight	0,9 kg
Shipping method	Parcel service
Net weight approx.	0,90 kg
Gross weight approx.	0,90 kg
Shipping weight	0,9 kg

## Services (optional)

Article number for DAkkS calibration	963-216
--------------------------------------	---------