

# KERN 304-036

# KERN

Weight set, E1, 1 g - 100 g, Knob, Stainless steel polished, in aluminium case



## Pictograms

### STANDARD



## Category

|                  |            |
|------------------|------------|
| Brand            | KERN       |
| Product category | Weight set |
| Product group    | OIML E1    |
| Product family   | 304        |

## Measuring System

|                             |             |
|-----------------------------|-------------|
| Test weight - nominal value | 1 g - 100 g |
| Weight OIML class           | E1          |

## Construction

|                    |                          |
|--------------------|--------------------------|
| Shape weight       | Knob                     |
| Dimension (W×D×H)  | 165×125×75 mm            |
| Case of weight set | in aluminium case        |
| Material           | Stainless steel polished |

## Packing & Shipping

|                              |                |
|------------------------------|----------------|
| Readability force [d] (N)    | 1 d            |
| Dimensions packaging (W×D×H) | 165×125×75 mm  |
| Net weight                   | 0,65 kg        |
| Shipping method              | Parcel service |
| Net weight approx.           | 0,65 kg        |
| Gross weight approx.         | 0,90 kg        |
| Shipping weight              | 0,9 kg         |

## Services (optional)

|  |         |
|--|---------|
| Article number for DAkkS calibration         | 963-216 |
| Article number for guaranteed plus-deviation | 921-216 |