

KERN 304-04

KERN

Weight set, E1, 1 g - 200 g, Knob, Stainless steel polished, in wooden case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E1
Product family	304

Measuring System

Test weight - nominal value	1 g - 200 g
Weight OIML class	E1

Construction

Shape weight	Knob
Dimension (W×D×H)	152×152×66 mm
Number of weights	11
Case of weight set	in wooden case
Material	Stainless steel polished

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	152×152×66 mm
Net weight	1,27 kg
Shipping method	Parcel service
Net weight approx.	1,4 kg
Gross weight approx.	1,4 kg
Shipping weight	1,27 kg

Services (optional)

Article number for DAkkS calibration	963-217
Article number for guaranteed plus-deviation	921-217