

# KERN 304-046

# KERN

Weight set, E1, 1 g - 200 g, Knob, Stainless steel polished, in aluminium case



## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML E1
Product family	304

## Measuring System

Test weight - nominal value	1 g - 200 g
Weight OIML class	E1

## Construction

Shape weight	Knob
Dimension (W×D×H)	190×130×80 mm
Case of weight set	in aluminium case
Material	Stainless steel polished

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	190×130×80 mm
Net weight	0,9 kg
Shipping method	Parcel service
Net weight approx.	0,90 kg
Gross weight approx.	1,2 kg
Shipping weight	1,1 kg

## Services (optional)

Article number for DAkkS calibration	963-217
Article number for guaranteed plus-deviation	921-217