

KERN 304-054

KERN

Weight set, E1, 1 g - 500 g, Knob, Stainless steel polished, in plastic case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E1
Product family	304

Measuring System

Test weight - nominal value	1 g - 500 g
Weight OIML class	E1

Construction

Shape weight	Knob
Dimension (W×D×H)	250×210×75 mm
Number of weights	12
Case of weight set	in plastic case
Material	Stainless steel polished

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	250×210×75 mm
Net weight	3,19 kg
Shipping method	Parcel service
Net weight approx.	3,2 kg
Gross weight approx.	3,2 kg
Shipping weight	3,19 kg

Services (optional)

Article number for DAkkS calibration	963-218
--------------------------------------	---------