

KERN 304-06

KERN

Weight set, E1, 1 g - 1 kg, Knob, Stainless steel polished, in wooden case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E1
Product family	304

Measuring System

Test weight - nominal value	1 g - 1 kg
Weight OIML class	E1

Construction

Shape weight	Knob
Dimension (W×D×H)	230×180×110 mm
Number of weights	13
Case of weight set	in wooden case
Material	Stainless steel polished

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	230×180×110 mm
Net weight	4,21 kg
Shipping method	Parcel service
Net weight approx.	4,4 kg
Gross weight approx.	4,4 kg
Shipping weight	4,21 kg

Services (optional)

Article number for DAkkS calibration	963-219
Article number for guaranteed plus-deviation	921-219