

KERN 304-076

KERN

Weight set, E1, 1 g - 2 kg, Knob, Stainless steel polished, in aluminium case



Pictograms

STANDARD



Category

| | |
|------------------|------------|
| Brand | KERN |
| Product category | Weight set |
| Product group | OIML E1 |
| Product family | 304 |

Measuring System

| | |
|-----------------------------|------------|
| Test weight - nominal value | 1 g - 2 kg |
| Weight OIML class | E1 |

Construction

| | |
|--------------------|--------------------------|
| Shape weight | Knob |
| Dimension (W×D×H) | 278×180×138 mm |
| Case of weight set | in aluminium case |
| Material | Stainless steel polished |

Packing & Shipping

| | |
|------------------------------|----------------|
| Readability force [d] (N) | 1 d |
| Dimensions packaging (W×D×H) | 278×180×138 mm |
| Net weight | 5,15 kg |
| Shipping method | Parcel service |
| Net weight approx. | 6 kg |
| Gross weight approx. | 6 kg |
| Shipping weight | 5,55 kg |

Services (optional)

| | |
|--|---------|
| Article number for DAkkS calibration | 963-220 |
| Article number for guaranteed plus-deviation | 921-220 |