

KERN 304-084

KERN

Weight set, E1, 1 g - 5 kg, Knob, Stainless steel polished, in plastic case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E1
Product family	304

Measuring System

Test weight - nominal value	1 g - 5 kg
Weight OIML class	E1

Construction

Shape weight	Knob
Dimension (W×D×H)	350×320×180 mm
Number of weights	16
Case of weight set	in plastic case
Material	Stainless steel polished

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	350×320×180 mm
Net weight	12,69 kg
Shipping method	Parcel service
Net weight approx.	13 kg
Gross weight approx.	13 kg
Shipping weight	12,69 kg

Services (optional)

Article number for DAkkS calibration	963-221
--------------------------------------	---------