

# KERN 307-09

# KERN

Individual weight, OIML Class E1, 500 g, Knob, Stainless steel polished



## Services (optional)

Article number for DAkkS calibration	963-239
Article number for guaranteed plus-deviation	921-239

## Pictograms

### STANDARD



### OPTION



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML E1
Product family	307

## Measuring System

Test weight - nominal value	500 g
Test weight - volume (cm <sup>3</sup> )	62,5
Test weight - volume uncertainty (cm <sup>3</sup> )	0,7813
Permissible tolerance according to OIML (+/-)	0,25 mg
Test weight - Density [Min] (OIML)	7934
Test weight - Density [Max] (OIML)	8067
Weight OIML class	E1

## Construction

Shape weight	Knob
Dimensions (Ø×H)	38×63 mm
Material	Stainless steel polished

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	38×38×63 mm
Net weight	0,5 kg
Shipping method	Parcel service
Net weight approx.	0,50 kg
Gross weight approx.	0,50 kg
Shipping weight	0,5 kg