

KERN 307-11

KERN

Individual weight, OIML Class E1, 1 kg, Knob, Stainless steel polished



Services (optional)

Article number for DAkkS calibration	963-241
Article number for guaranteed plus-deviation	921-241

Pictograms

STANDARD



OPTION



Category

Brand	KERN
Product category	Individual weight
Product group	OIML E1
Product family	307

Measuring System

Test weight - nominal value	1 kg
Test weight - volume (cm ³)	125
Test weight - volume uncertainty (cm ³)	1,5625
Permissible tolerance according to OIML (+/-)	0,5 mg
Test weight - Density [Min] (OIML)	7934
Test weight - Density [Max] (OIML)	8067
Weight OIML class	E1

Construction

Shape weight	Knob
Dimensions (Ø×H)	48×80 mm
Material	Stainless steel polished

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	48×48×80 mm
Net weight	1 kg
Shipping method	Parcel service
Net weight approx.	1,00 kg
Gross weight approx.	1,00 kg
Shipping weight	1 kg