

KERN 307-14

Individual weight, OIML Class E1, 10 kg, Knob, Stainless steel polished



Services (optional)

Article number for DAkkS calibration	963-244
Article number for guaranteed plus-deviation	921-244

Pictograms

STANDARD



OPTION



Category

Brand	KERN
Product category	Individual weight
Product group	OIML E1
Product family	307

Measuring System

Test weight - nominal value	10 kg
Test weight - volume (cm ³)	1250
Test weight - volume uncertainty (cm ³)	15,625
Permissible tolerance according to OIML (+/-)	5 mg
Test weight - Density [Min] (OIML)	7934
Test weight - Density [Max] (OIML)	8067
Weight OIML class	E1

Construction

Shape weight	Knob
Dimensions (Ø×H)	100×185 mm
Material	Stainless steel polished

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	100×100×185 mm
Net weight	10 kg
Shipping method	Parcel service
Net weight approx.	10 kg
Gross weight approx.	10 kg
Shipping weight	10 kg