

KERN 307-16

Individual weight, OIML Class E1, 50 kg, Knob, Stainless steel polished



Services (optional)

Article number for DAkkS calibration	963-246
Article number for guaranteed plus-deviation	921-246

Pictograms

STANDARD



OPTION



Category

Brand	KERN
Product category	Individual weight
Product group	OIML E1
Product family	307

Measuring System

Test weight - nominal value	50 kg
Test weight - volume (cm ³)	6250
Test weight - volume uncertainty (cm ³)	78,125
Permissible tolerance according to OIML (+/-)	25 mg
Test weight - Density [Min] (OIML)	7934
Test weight - Density [Max] (OIML)	8067
Weight OIML class	E1

Construction

Shape weight	Knob
Dimensions (Ø×H)	178×298 mm
Material	Stainless steel polished

Packing & Shipping

Readability force [d] (N)	2 d
Dimensions packaging (W×D×H)	178×178×298 mm
Net weight	50 kg
Shipping method	Freight
Net weight approx.	50 kg
Gross weight approx.	50 kg
Shipping weight	50 kg