

KERN 308-37

KERN

Milligram weights, E1, 100 mg, Wire, stainless steel polished (OIML)



Services (optional)

Article number for DAkkS calibration 962-257

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML E1
Product family	308

Measuring System

Test weight - nominal value	100 mg
Test weight - volume (cm ³)	0,0125
Test weight - volume uncertainty (cm ³)	0,0003
Permissible tolerance according to OIML (+/-)	0,005 mg
Test weight - Density [Min] (OIML)	4400
Test weight - Density [Max] (OIML)	100000
Weight OIML class	E1

Construction

Shape weight	Wire
Dimensions weight (W×D×H)	9×9×3 mm
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	9×9×3 mm
Net weight	0,0001 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg