

KERN 308-42

KERN

Set of milligram weights, E1, 1 mg - 500 mg, Wire, Stainless steel polished, in plastic case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E1
Product family	308

Measuring System

Test weight - nominal value	1 mg - 500 mg
Weight OIML class	E1

Construction

Shape weight	Wire
Dimension (W×D×H)	136×53×22 mm
Number of weights	13
Case of weight set	in plastic case
Material	Stainless steel polished

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	136×53×22 mm
Net weight	0,09 kg
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,10 kg
Shipping weight	0,1 kg

Services (optional)

Article number for DAkkS calibration	962-250
--------------------------------------	---------