

KERN 312-026

KERN

Weight set, E2, 1 g - 50 g, Compact form, stainless steel polished (OIML), in aluminium case



Pictograms

STANDARD



Category

| | |
|------------------|------------|
| Brand | KERN |
| Product category | Weight set |
| Product group | OIML E2 |
| Product family | 312 |

Measuring System

| | |
|-----------------------------|------------|
| Test weight - nominal value | 1 g - 50 g |
| Weight OIML class | E2 |

Construction

| | |
|--------------------|---------------------------------|
| Shape weight | Compact form |
| Dimension (W×D×H) | 165×120×75 mm |
| Case of weight set | in aluminium case |
| Material | stainless steel polished (OIML) |

Packing & Shipping

| | |
|------------------------------|----------------|
| Readability force [d] (N) | 1 d |
| Dimensions packaging (W×D×H) | 165×120×75 mm |
| Net weight | 0,54 kg |
| Shipping method | Parcel service |
| Net weight approx. | 0,55 kg |
| Gross weight approx. | 0,80 kg |
| Shipping weight | 0,8 kg |

Services (optional)

| | |
|--|---------|
| Article number for DAkkS calibration | 962-315 |
| Article number for guaranteed plus-deviation | 921-315 |