

KERN 312-034

KERN

Weight set, E2, 1 g - 100 g, Compact form, stainless steel polished (OIML), in plastic case



Article number for guaranteed
plus-deviation

921-316

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E2
Product family	312

Measuring System

Test weight - nominal value	1 g - 100 g
OIML-compliant	✓
Weight OIML class	E2

Construction

Shape weight	Compact form
Dimension (W×D×H)	250×210×75 mm
Number of weights	9
Case of weight set	in plastic case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	250×210×75 mm
Net weight	0,6 kg
Shipping method	Parcel service
Net weight approx.	0,60 kg
Gross weight approx.	0,60 kg
Shipping weight	0,8 kg

Services (optional)

Article number for DAKkS calibration	962-316
---	---------