

# KERN 312-034

# KERN

Weight set, E2, 1 g - 100 g, Compact form, stainless steel polished (OIML), in plastic case



## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML E2
Product family	312

## Measuring System

Test weight - nominal value	1 g - 100 g
Weight OIML class	E2

## Construction

Shape weight	Compact form
Dimension (W×D×H)	250×210×75 mm
Number of weights	9
Case of weight set	in plastic case
Material	stainless steel polished (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	250×210×75 mm
Net weight	0,6 kg
Shipping method	Parcel service
Net weight approx.	0,60 kg
Gross weight approx.	0,60 kg
Shipping weight	0,8 kg

## Services (optional)

Article number for DAkkS calibration	962-316
Article number for guaranteed plus-deviation	921-316