

KERN 312-044

KERN

Weight set, E2, 1 g - 200 g, Compact form, stainless steel polished (OIML), in plastic case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E2
Product family	312

Measuring System

Test weight - nominal value	1 g - 200 g
Weight OIML class	E2

Construction

Shape weight	Compact form
Dimension (W×D×H)	250×210×75 mm
Number of weights	11
Case of weight set	in plastic case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	250×210×75 mm
Net weight	0,94 kg
Shipping method	Parcel service
Net weight approx.	0,95 kg
Gross weight approx.	0,95 kg
Shipping weight	0,94 kg

Services (optional)

Article number for DAkkS calibration	962-317
Article number for guaranteed plus-deviation	921-317