

KERN 312-066

KERN

Weight set, E2, 1 g - 1 kg, Compact form, stainless steel polished (OIML), in aluminium case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E2
Product family	312

Measuring System

Test weight - nominal value	1 g - 1 kg
Weight OIML class	E2

Construction

Shape weight	Compact form
Dimension (W×D×H)	230×140×110 mm
Case of weight set	in aluminium case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	230×140×110 mm
Net weight	2,7 kg
Shipping method	Parcel service
Net weight approx.	2,8 kg
Gross weight approx.	3,6 kg
Shipping weight	3,5 kg

Services (optional)

Article number for DAkkS calibration	962-319
Article number for guaranteed plus-deviation	921-319