

KERN 313-03

KERN

Weight set, E2, 1 mg - 100 g, Knob, stainless steel polished (OIML), in wooden case



Article number for guaranteed plus-deviation 921-302

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E2
Product family	313

Measuring System

Test weight - nominal value	1 mg - 100 g
Weight OIML class	E2

Construction

Shape weight	Knob
Dimension (W×D×H)	152×152×66 mm
Number of weights	21
Case of weight set	in wooden case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	152×152×66 mm
Net weight	0,9 kg
Shipping method	Parcel service
Net weight approx.	0,90 kg
Gross weight approx.	0,90 kg
Shipping weight	0,9 kg

Services (optional)

Article number for conformity assessment	952-302
Article number for DAkkS calibration	962-302