

# KERN 313-054

# KERN

Weight set, E2, 1 mg - 500 g, Knob, stainless steel polished (OIML), in plastic case



Article number for guaranteed plus-deviation 921-304

Pictograms

STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML E2
Product family	313

## Measuring System

Test weight - nominal value	1 mg - 500 g
Weight OIML class	E2

## Construction

Shape weight	Knob
Dimension (W×D×H)	250×210×75 mm
Number of weights	24
Case of weight set	in plastic case
Material	stainless steel polished (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	250×210×75 mm
Net weight	1,44 kg
Shipping method	Parcel service
Net weight approx.	1,6 kg
Gross weight approx.	1,6 kg
Shipping weight	1,44 kg

## Services (optional)

Article number for conformity assessment	952-304
Article number for DAkkS calibration	962-304