

KERN 313-06

KERN

Weight set, E2, 1 mg - 1 kg, Knob, stainless steel polished (OIML), in wooden case



Article number for guaranteed plus-deviation 921-305

Pictograms

STANDARD



Category

| | |
|------------------|------------|
| Brand | KERN |
| Product category | Weight set |
| Product group | OIML E2 |
| Product family | 313 |

Measuring System

| | |
|-----------------------------|-------------|
| Test weight - nominal value | 1 mg - 1 kg |
| Weight OIML class | E2 |

Construction

| | |
|--------------------|---------------------------------|
| Shape weight | Knob |
| Dimension (WxDxH) | 230x180x110 mm |
| Number of weights | 25 |
| Case of weight set | in wooden case |
| Material | stainless steel polished (OIML) |

Packing & Shipping

| | |
|------------------------------|----------------|
| Readability force [d] (N) | 1 d |
| Dimensions packaging (WxDxH) | 230x180x110 mm |
| Net weight | 4,21 kg |
| Shipping method | Parcel service |
| Net weight approx. | 4,4 kg |
| Gross weight approx. | 4,4 kg |
| Shipping weight | 4,21 kg |

Services (optional)

| | |
|--|---------|
| Article number for conformity assessment | 952-305 |
| Article number for DAkkS calibration | 962-305 |