

# KERN 313-060-100

Wooden box for weight sets (Cylindrical/compact, 1 mg - 1 kg)



## Category

|                  |               |
|------------------|---------------|
| Brand            | KERN          |
| Product category | Cases & Boxes |
| Product group    | Wooden box    |

## Measuring System

|   |                     |
|---|---------------------|
| Case suitable for OIML classes          | E1+E2+F1            |
| Case suitable for nominal weight values | 1 mg - 1 kg         |
| Suitable weight shape for case          | Cylindrical/compact |

## Construction

|                   |                |
|-------------------|----------------|
| Dimension (W×D×H) | 230×180×110 mm |
| Material          | beech          |

## Packing & Shipping

|                              |                |
|------------------------------|----------------|
| Readability force [d] (N)    | 1 d            |
| Dimensions packaging (W×D×H) | 230×180×110 mm |
| Net weight                   | 2,2 kg         |
| Shipping method              | Parcel service |
| Net weight approx.           | 2,2 kg         |
| Gross weight approx.         | 2,6 kg         |
| Shipping weight              | 2,5 kg         |

## Pictograms

STANDARD

