

KERN 313-060-600

KERN

Aluminium protective box for weight sets (Cylindrical/compact, 1 mg - 1 kg)



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Aluminium protective box

Measuring System

Case suitable for OIML classes	E1 - M1
Case suitable for nominal weight values	1 mg - 1 kg
Suitable weight shape for case	Cylindrical/compact

Construction

Dimension (W×D×H)	230×140×110 mm
Material	aluminium

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	230×140×110 mm
Net weight	0,8 kg
Shipping method	Parcel service
Net weight approx.	0,80 kg
Gross weight approx.	0,80 kg
Shipping weight	0,8 kg

Pictograms

STANDARD

