

KERN 313-064

KERN

Weight set, E2, 1 mg - 1 kg, Knob, stainless steel polished (OIML), in plastic case



Article number for guaranteed plus-deviation 921-305

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E2
Product family	313

Measuring System

Test weight - nominal value	1 mg - 1 kg
Weight OIML class	E2

Construction

Shape weight	Knob
Dimension (W×D×H)	350×320×180 mm
Number of weights	25
Case of weight set	in plastic case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	350×320×180 mm
Net weight	3,7 kg
Shipping method	Parcel service
Net weight approx.	3,8 kg
Gross weight approx.	3,8 kg
Shipping weight	4 kg

Services (optional)

Article number for conformity assessment	952-305
Article number for DAkkS calibration	962-305