

KERN 313-07

KERN

Weight set, E2, 1 mg - 2 kg, Knob, stainless steel polished (OIML), in wooden case



Article number for guaranteed plus-deviation 921-306

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E2
Product family	313

Measuring System

Test weight - nominal value	1 mg - 2 kg
Weight OIML class	E2

Construction

Shape weight	Knob
Dimension (WxDxH)	305x190x140 mm
Number of weights	27
Case of weight set	in wooden case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (WxDxH)	305x190x140 mm
Net weight	9,01 kg
Shipping method	Parcel service
Net weight approx.	10 kg
Gross weight approx.	10 kg
Shipping weight	9,01 kg

Services (optional)

Article number for conformity assessment	952-306
Article number for DAkkS calibration	962-306