

KERN 313-074

KERN

Weight set, E2, 1 mg - 2 kg, Knob, stainless steel polished (OIML), in plastic case



Article number for guaranteed plus-deviation 921-306

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E2
Product family	313

Measuring System

Test weight - nominal value	1 mg - 2 kg
Weight OIML class	E2

Construction

Shape weight	Knob
Dimension (WxDxH)	350x320x180 mm
Number of weights	27
Case of weight set	in plastic case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (WxDxH)	350x320x180 mm
Net weight	7,69 kg
Shipping method	Parcel service
Net weight approx.	8 kg
Gross weight approx.	8 kg
Shipping weight	7,69 kg

Services (optional)

Article number for conformity assessment	952-306
Article number for DAkkS calibration	962-306