

KERN 313-086

KERN

Weight set, E2, 1 mg - 5 kg, Knob, stainless steel polished (OIML), in aluminium case



Article number for guaranteed plus-deviation 921-307

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E2
Product family	313

Measuring System

Test weight - nominal value	1 mg - 5 kg
Weight OIML class	E2

Construction

Shape weight	Knob
Dimension (W×D×H)	460×260×200 mm
Case of weight set	in aluminium case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	460×260×200 mm
Net weight	18 kg
Shipping method	Parcel service
Net weight approx.	18 kg
Gross weight approx.	24 kg
Shipping weight	23 kg

Services (optional)

Article number for conformity assessment	952-307
Article number for DAKkS calibration	962-307