

KERN 313-096

KERN

Weight set, E2, 1 mg - 10 kg, Knob, stainless steel polished (OIML), in aluminium case



Article number for guaranteed plus-deviation 921-308

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E2
Product family	313

Measuring System

Test weight - nominal value	1 mg - 10 kg
Weight OIML class	E2

Construction

Shape weight	Knob
Dimension (W×D×H)	610×260×235 mm
Case of weight set	in aluminium case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	610×260×235 mm
Net weight	23,99 kg
Shipping method	Parcel service
Net weight approx.	24 kg
Gross weight approx.	28 kg
Shipping weight	28 kg

Services (optional)

Article number for conformity assessment	952-308
Article number for DAKkS calibration	962-308