

# KERN 314-046

# KERN

Weight set, E2, 1 g - 200 g, Knob, stainless steel polished (OIML), in aluminium case



Article number for guaranteed plus-deviation 921-317

## Pictograms

## STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML E2
Product family	314

## Measuring System

Test weight - nominal value	1 g - 200 g
Weight OIML class	E2

## Construction

Shape weight	Knob
Dimension (W×D×H)	190×130×80 mm
Case of weight set	in aluminium case
Material	stainless steel polished (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	190×130×80 mm
Net weight	0,9 kg
Shipping method	Parcel service
Net weight approx.	0,90 kg
Gross weight approx.	1,2 kg
Shipping weight	1,1 kg

## Services (optional)

Article number for conformity assessment	952-317
Article number for DAkkS calibration	962-317