

KERN 314-054

KERN

Weight set, E2, 1 g - 500 g, Knob, stainless steel polished (OIML), in plastic case



Article number for guaranteed plus-deviation 921-318

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E2
Product family	314

Measuring System

Test weight - nominal value	1 g - 500 g
Weight OIML class	E2

Construction

Shape weight	Knob
Dimension (WxDxH)	250x210x75 mm
Number of weights	12
Case of weight set	in plastic case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (WxDxH)	250x210x75 mm
Net weight	1,44 kg
Shipping method	Parcel service
Net weight approx.	1,6 kg
Gross weight approx.	1,6 kg
Shipping weight	1,44 kg

Services (optional)

Article number for conformity assessment	952-318
Article number for DAkkS calibration	962-318