

KERN 314-066

KERN

Weight set, E2, 1 g - 1 kg, Knob, stainless steel polished (OIML), in aluminium case



Article number for guaranteed plus-deviation 921-319

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E2
Product family	314

Measuring System

Test weight - nominal value	1 g - 1 kg
Weight OIML class	E2

Construction

Shape weight	Knob
Dimension (W×D×H)	230×140×110 mm
Case of weight set	in aluminium case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	230×140×110 mm
Net weight	2,7 kg
Shipping method	Parcel service
Net weight approx.	2,8 kg
Gross weight approx.	3,6 kg
Shipping weight	3,5 kg

Services (optional)

Article number for conformity assessment	952-319
Article number for DAKkS calibration	962-319