

KERN 314-08

KERN

Weight set, E2, 1 g - 5 kg, Knob, stainless steel polished (OIML), in wooden case



Article number for guaranteed plus-deviation 921-321

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E2
Product family	314

Measuring System

Test weight - nominal value	1 g - 5 kg
Weight OIML class	E2

Construction

Shape weight	Knob
Dimension (W×D×H)	420×200×200 mm
Number of weights	16
Case of weight set	in wooden case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	420×200×200 mm
Net weight	17,81 kg
Shipping method	Parcel service
Net weight approx.	18 kg
Gross weight approx.	18 kg
Shipping weight	17,81 kg

Services (optional)

Article number for conformity assessment	952-321
Article number for DAkkS calibration	962-321