

# KERN 314-09

# KERN

Weight set, E2, 1 g - 10 kg, Knob, stainless steel polished (OIML), in wooden case



Article number for guaranteed plus-deviation 921-322

Pictograms

STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML E2
Product family	314

## Measuring System

Test weight - nominal value	1 g - 10 kg
Weight OIML class	E2

## Construction

Shape weight	Knob
Dimension (W×D×H)	540×200×240 mm
Number of weights	17
Case of weight set	in wooden case
Material	stainless steel polished (OIML)

## Packing & Shipping

Readability force [d] (N)	2 d
Dimensions packaging (W×D×H)	540×200×240 mm
Net weight	30,61 kg
Shipping method	Freight
Net weight approx.	32 kg
Gross weight approx.	34 kg
Shipping weight	32,5 kg

## Services (optional)

Article number for conformity assessment	952-322
Article number for DAkkS calibration	962-322