

# KERN 315-040-100



Wooden box for individual weight sets



## Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

## Measuring System

Case suitable for OIML classes	E1+E2+F1
Case suitable for nominal weight values	universal bis 200 g
Suitable weight shape for case	Cylindrical/compact

## Construction

Dimension (W×D×H)	152×152×66 mm
Material	beech

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	150×140×65 mm
Net weight	0,68 kg
Shipping method	Parcel service
Net weight approx.	0,70 kg
Gross weight approx.	0,70 kg
Shipping weight	0,68 kg

## Pictograms

STANDARD

