

# KERN 315-109-205

Safety Set 20 g (F1); 500 g (E2)

# KERN



## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	Safety Set
Product family	315

## Measuring System

Test weight - nominal value	20 g & 500 g
Weight OIML class	E2 F1

## Construction

Number of weights	2
Case of weight set	in aluminium case
Material	stainless steel

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	255×135×140 mm
Net weight	1,2 kg
Shipping method	Parcel service
Net weight approx.	1,2 kg
Gross weight approx.	3,0 kg
Shipping weight	3 kg

## Services (optional)

Article number for DAkkS calibration	962-300-400
--------------------------------------	-------------