

KERN 315-111-206

Safety Set 50 g (F1); 1000 g (E2)

KERN



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	Safety Set
Product family	315

Measuring System

Test weight - nominal value	50 g & 1 kg
Weight OIML class	E2 F1

Construction

Number of weights	2
Case of weight set	in aluminium case
Material	stainless steel

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	255×135×140 mm
Net weight	1,2 kg
Shipping method	Parcel service
Net weight approx.	1,2 kg
Gross weight approx.	3,0 kg
Shipping weight	3 kg

Services (optional)

Article number for DAkkS calibration	962-300-400
--------------------------------------	-------------