

KERN 315-243

KERN

Tweezers to be able to safely grip small test weights



Category

Brand	KERN
Product category	Test weights accessories
Product group	Tweezer

Construction

Dimension (W×D×H)	103×5×8 mm
-------------------	------------

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	105×10×10 mm
Net weight	0,02 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

Pictograms

STANDARD

