

KERN 315-275

Gewicht handvat, plastic gecoat

KERN



Categorie

Merk	KERN
Productcategorie	Controlegewichtenaccessoires
Productgroep	Gewicht handvat
Productfamilie	315

Verpakking & verzending

Afreesbaarheid kracht [d] (N)	1 d
Afmetingen verpakking (B×D×H)	90×160×95 mm
Nettogewicht	0,35 kg
Verzendmethode	Pakketdienst
Nettogewicht ca.	0,35 kg
Brutogewicht ca.	0,45 kg
Verzendgewicht	0,45 kg

Pictogrammen

STANDARD

