

KERN 315-309-205

Safety Set 500 g (F2); 20 g (F1)

KERN



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	Safety Set
Product family	315

Measuring System

Test weight - nominal value	20 g & 500 g
Weight OIML class	F1 F2

Construction

Number of weights	2
Case of weight set	in aluminium case
Material	stainless steel

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	255×135×140 mm
Net weight	1,2 kg
Shipping method	Parcel service
Net weight approx.	1,2 kg
Gross weight approx.	3,0 kg
Shipping weight	3 kg

Services (optional)

Article number for DAkkS calibration	962-400-400
Article number for guaranteed plus-deviation	921-439